

Demonstration Video: Injecting with the SureClick® Autoinjector

VO:

LORIE: Oh! Hey there, hand tight - the SureClick® demo video will begin after this brief safety message.

INDICATIONS:

Prescription Enbrel® (etanercept) is taken or given by injection.

The video you are about to see is intended to provide you with an overview and basic information about the ENBREL SureClick® Autoinjector. This video is not a comprehensive review of the ENBREL SureClick® Autoinjector and does not replace the Instructions for Use. Please review with your healthcare provider the Patient Instructions for Use included with ENBREL SureClick® Autoinjector. Be sure to talk to your healthcare provider about any questions or concerns you may have about the ENBREL SureClick® Autoinjector.

The Important Safety Information, including Boxed Warning for ENBREL, will be discussed later in this video. ENBREL is indicated for reducing signs and symptoms, keeping joint damage from getting worse, and improving physical function in patients with moderately to severely active rheumatoid arthritis. ENBREL can be taken with methotrexate or used alone.

ENBREL is indicated for reducing signs and symptoms of moderately to severely active polyarticular juvenile idiopathic arthritis (pJIA) in children 2 years of age or older.

ENBREL is indicated for reducing signs and symptoms, keeping joint damage from getting worse, and improving physical function in adult patients with psoriatic arthritis. ENBREL can be used with or without methotrexate.

ENBREL is indicated for active juvenile psoriatic arthritis (JPsA) in children 2 years of age or older.

ENBREL is indicated for reducing signs and symptoms in patients with active ankylosing spondylitis.

ENBREL is indicated for chronic moderate to severe plaque psoriasis (PsO) in children 4 years of age or older and adults who may benefit from taking injections or pills (systemic therapy) or phototherapy (ultraviolet light).

Individual results may vary. Talk to your doctor about the possible risks and benefits of ENBREL.

Please see Important Safety Considerations later in this video and the Important Safety Information, including Boxed Warning, on this webpage. Please also see the full Prescribing Information and Medication Guide available on this website.

MARY: Coming!

LORIE: Mary?

MARY: Hi!

LORIE: Hi, I'm Lorie, your ENBREL Nurse Partner.

MARY: Hi, Lorie. Come on in! Thanks for coming over.

LORIE: Absolutely, my pleasure.

LORIE: So as you know, I'm here to help you as you get started with the SureClick® Autoinjector.

MARY: Great.

LORIE: Think of today as a demonstration that can help supplement the training that you received in your physician's office. So, I'll go over things in detail. But before I do, your SureClick® came with Instructions for Use.

MARY: I have one, yup. I reviewed it with the nurse when we received training in the office and read it just before you came over. I'm really glad you are here, because this is the first time I'm injecting myself at home.

LORIE: Oh, you're welcome and that's great. You should always read the instructions for use before injecting and also refer to it if you have any questions.

MARY: Great.

LORIE: Alright, so are you ready to get started?

MARY: Ready.

LORIE: Now, before you do the actual injection, I want to show you some features of the SureClick® Autoinjector.

LORIE: So, the first thing is understanding which end is which. So, see this purple button? That's the top, and that's the start button. Your thumb goes here. This white cap covers the needle. Do not remove the white cap until you are actually ready to inject. And now, once you remove this white cap... you'll see a green safety guard. This is where the needle is tucked away inside. And that's the end you want to place against your skin.

MARY: I don't see the needle.

LORIE: Exactly!

MARY: Oh.

LORIE: The SureClick® is designed so you don't have to see the needle. So, keep in mind, once you remove this white cap, you want to inject within 5 minutes to prevent the medicine from drying out. Now there's one more part I want you to look at. See this window? That clear liquid is the ENBREL medicine. Always make sure that it's clear, has no clumps or flakes floating in it. If the medicine inside looks cloudy, discolored or your SureClick® has expired, don't use it. You'll want to contact your doctor right away. Or you can call 1-888-4ENBREL to get a replacement. So, where's your ENBREL SureClick® now?

MARY: In my fridge. I put it in there right after I got it.

LORIE: Perfect. Let's go get one of the SureClick®s out of the box.

MARY: Okay, follow me. I need to warn you, my fridge is a bit messy. Don't judge me.

LORIE: Don't worry! You should see my fridge!

LORIE: I see you keep your SureClick® on the top shelf. That's perfect, you want to keep it out of reach of children.

MARY: Got it.

LORIE: Couple quick things to keep in mind. Never freeze it or store it in extreme heat or cold. Don't shake it. And always store it in its original carton, like you have. That protects it from light damage. And remember, leave your unused medication in the refrigerator. Now take your SureClick® out of the fridge then wait. You want this to get to room temperature. Injecting at room temperature can make for a more comfortable injection. So, let it sit out for about 30 minutes. While you wait, you can stream a show. Take a walk. Just make sure it's out of reach of small children and pets.

MARY: Ah, I know just the place. Oh, question - I'm going on vacation soon. I'm going to be staying at a hotel for about a week, and I'm not sure if they have a refrigerator or anything. Can I travel with ENBREL?

LORIE: Absolutely. ENBREL can actually be at room temperature for up to 30 days. So you can actually bring it with you in your carry-on, and keep taking it right on schedule.

MARY: That's perfect.

LORIE: But for today, let's just do 30 minutes.

MARY: Okay.

LORIE: Oh, one more thing. Once ENBREL comes to room temperature, do not put it back in the refrigerator.

MARY: Good to know. Wow. That was a fast 30 minutes.

LORIE: Well, first things first... let's wash your hands thoroughly.

MARY: I'm all set.

LORIE: Alright, let's go get your SureClick® down.

MARY: Okay.

LORIE: Now, this is everything you'll need. Cotton balls, alcohol wipes, a bandage, and a sharps container for when you're done.

MARY: Great.

LORIE: So, today you'll be injecting into your thigh.

MARY: Okay.

LORIE: But you can see other injection areas in your Instructions for Use.

LORIE: Now, let's use an alcohol wipe to clean a 2-inch circle around where you're going to inject.

MARY: Okay.

LORIE: And then we'll let it dry. Okay, now hold the SureClick® so you can easily see this window when you're injecting.

MARY: Okay.

LORIE: And you'll want to remove the white cap. And remember don't remove the white cap more than 5 minutes prior to injecting to avoid the medicine drying out. So, now you'll want to stretch your skin slightly to create a firm surface. There you go, perfect. Hold the SureClick® at a 90-degree angle on your skin. An easy way to remember these next steps is something I call "push, press, count". Now push the green safety guard firmly against your skin - remembering 90 degrees. Now this is important - Don't touch the purple button until you're actually ready to inject.

MARY: I think I'm ready.

LORIE: Okay! And remember keep pushing the SureClick® against your skin. It takes about 15 seconds to get the full dose.

MARY: Okay.

LORIE: So, we're going slowly to count slowly to 15, just to make sure. Okay, now press the purple button all the way down and release the button. And now count...

MARY: One, two, three, four, five...

LORIE: And don't forget to keep an eye on that window, too.

MARY: ...twelve, thirteen, fourteen, fifteen. The window is yellow now.

LORIE: Great! That means the injection is complete. So that's an important reason why you always want to check that window before you remove it from the skin. If you remove the

SureClick® too early, you may not receive your full dose of ENBREL. So, call your doctor if you think you didn't receive the full dose. And remember: Push. Press. Count.

MARY: Push. Press. Count. Thanks Lorie. It was helpful having you here today.

LORIE: Oh, that's great! So, there's a bandage there if you need it.

MARY: Great.

LORIE: And next, you're going to want to dispose of your empty SureClick® safely. There's no medicine in it anymore, but it does have a needle, so you can't just throw them in your household trash or recycling. You want to put the used SureClick® autoinjector in an FDA-cleared sharps disposal container right away after use. Drop it in, seal it up and make sure your kids and pets can't get to it.

MARY: That's it?

LORIE: That's it.

MARY: Thanks so much for coming. I'm really much more comfortable now that we went over this.

LORIE: Of course, it's my pleasure! And remember to refer to your Instructions for Use if you have any questions. And you can also call Enbrel Support. We're happy to help!

MARY: Awesome. Bye now!

LORIE: Bye!

IMPORTANT SAFETY INFORMATION

What is the most important information I should know about ENBREL?

ENBREL is a medicine that affects your immune system. ENBREL can lower the ability of your immune system to fight infections. Serious infections have happened in patients taking ENBREL. These infections include tuberculosis (TB) and infections caused by viruses, fungi, or bacteria that have spread throughout the body. Some patients have died from these infections. Your healthcare provider should test you for TB before you take ENBREL and monitor you closely for TB before, during, and after ENBREL treatment, even if you have tested negative for TB.

There have been some cases of unusual cancers, some resulting in death, reported in children and teenagers who started using tumor necrosis factor (TNF) blockers before 18 years of age. Also, for children, teenagers, and adults taking TNF blockers, including ENBREL, the chances of getting lymphoma or other cancers may increase. Patients with RA may be more likely to get lymphoma.

Before starting ENBREL, tell your healthcare provider if you:

- Have any existing medical conditions
- Are taking any medicines, including herbals
- Think you have, are being treated for, have signs of, or are prone to infection. You should not start taking ENBREL if you have any kind of infection, unless your healthcare provider says it is okay
- Have any open cuts or sores
- Have diabetes, HIV, or a weak immune system
- Have TB or have been in close contact with someone who has had TB
- Were born in, lived in, or traveled to countries where there is more risk for getting TB. Ask your healthcare provider if you are not sure
- Live, have lived in, or traveled to certain parts of the country (such as, the Ohio and Mississippi River valleys, or the Southwest) where there is a greater risk for certain kinds of fungal infections, such as histoplasmosis. These infections may develop or become more severe if you take ENBREL. If you don't know if these infections are common in the areas you've been to, ask your healthcare provider
- Have or have had hepatitis B
- Have or have had heart failure
- Develop symptoms such as persistent fever, bruising, bleeding, or paleness while taking ENBREL
- Use the medicine Kineret (anakinra), Orencia (abatacept), or Cytosan (cyclophosphamide)
- Are taking anti-diabetic medicines
- Have, have had, or develop a serious nervous disorder, seizures, any numbness or tingling, or a disease that affects your nervous system such as multiple sclerosis or Guillain-Barré syndrome
- Are scheduled to have surgery
- Have recently received or are scheduled for any vaccines. All vaccines should be brought up-to-date before starting ENBREL. Patients taking ENBREL should not receive live vaccines
- Are pregnant, planning to become pregnant, or breastfeeding
- Have been around someone with chicken pox

What are the possible side effects of ENBREL?

ENBREL can cause serious side effects including: New **infections** or worsening of infections you already have; **hepatitis B** can become active if you already have had it; **nervous system problems**, such as multiple sclerosis, seizures, or inflammation of the nerves of the eyes; **blood problems** (some fatal); new or worsening **heart failure**; new or worsening psoriasis; **allergic reactions**; **autoimmune reactions**, including a lupus-like syndrome and autoimmune hepatitis.

Common side effects include: Injection site reactions and upper respiratory infections (sinus infections).

In general, side effects in children were similar in frequency and type as those seen in adult patients. The types of infections reported were generally mild and similar to those usually seen in children.

These are not all the side effects with ENBREL. Tell your healthcare provider about any side effect that bothers you or does not go away.

If you have any questions about this information, be sure to discuss them with your healthcare provider. You are encouraged to report negative side effects of prescription drugs to the FDA. Visit www.fda.gov/medwatch, or call 1-800-FDA-1088.

Please see [Prescribing Information](#) and [Medication Guide](#) available on this website.



Manufactured by Immunex Corporation
Marketed by Amgen, Inc.

© 2022, 2024 Amgen, Inc. Thousand Oaks, CA 91320. All rights reserved. USA-916-83980